



ORTHOPEDICS

PROSTHETICS

SKIN CARE





OUR STORY

Founded in 1989, Silipos has established itself as the global leader in gel technologies across a wide range of industries including healthcare, consumer goods, and health and beauty aids. Our products are manufactured in Niagara Falls, NY and sold to over 7,000 customers in 90 countries worldwide.

Silipos is committed to expanding and improving our range of products to provide our customers with new and innovative gel solutions. All gel development including the incorporation of actives/additives is conducted in-house, as well as the base gel manufacturing and conversion into finished products. As an ISO 9001:2008 and ISO 13485:2003 certified company, we pride ourselves on product excellence and quality. Additionally, many of our products are CE marked for distribution in worldwide markets. Within the healthcare industry, Silipos is the most trusted source for medical solutions among healthcare professionals. We have built our reputation on innovation, quality, and outstanding customer service.

SOME OF THE MEDICAL SPECIALTIES TO WHICH WE PROVIDE SERVICE ARE:

- Podiatry/Footcare
- Orthopedics
- Prosthetics
- Physical Therapy
- Occupational Therapy
- Sports Medicine
- Scar Management
- Plastic Surgery
- Dermatology
- Health & Beauty/Skin Care
- Rheumatology
- Home Healthcare



OUR GELS

Silipos is the leading, vertically integrated manufacturer of patented polymer gel technologies in the world. With patented technology and extensive expertise in the manufacturing, processing, and conversion of gels, many leading global corporations rely on Silipos for their product development needs.

The unique “visco-elastic” properties of Silipos’ medical grade mineral oil gels allow them to move synergistically with the skin helping to prevent abrasion and discomfort. Additionally, Silipos gels are non-toxic, latex free and hypoallergenic.

Silipos produces over 200 products for a variety of product applications and uses. Silipos gels cushion and protect against impact, pressure, friction, and shear forces while simultaneously releasing moisture-rich medical grade mineral oil that gently hydrates, softens, soothes, and conditions the skin. Silipos gels have the unique ability to be formulated with a wide range of additives such as moisturizers, fragrances, vitamins, and minerals for a number of proven indications such as anti-microbial/anti-bacterial, antioxidants, anti-aging, and scar reduction.

Silipos gels are washable and reusable, and with proper care will continue to release skin beneficial ingredients throughout the life of the product.

Orthopedics


Antibacterial	6
Digital Care	9
Digital Care Kits	13
Bunion Care	14
Metatarsal	15
Sleeves & Wraps	18
Post Wound/Burn Care Kit	22
Heel Cups & Insoles	23
Gel Socks	26
Gel Pads & Sheeting	27
TheraStep®	30

Prosthetics

Sheaths	34
Liners	36
Sleeves	44
Pads & Accessories	46

Skin Care

Gel Care® Advanced	50
Scar Management	52
Gel Gloves & Socks	53
GeLuscious®	54





Orthopedics

Silipos' orthopedic and foot care products are made with medical grade mineral oil gels that easily conform to the body ensuring a comfortable fit that cushions and protects. When applied directly to the skin, our gels have the added benefit of skin moisturizing and conditioning.

Our orthopedic and foot care products offer easy to use solutions for a variety of ailments ranging from common foot discomforts to post traumatic injuries. Whatever your needs, Silipos has a product that can help improve the quality of life for your patients.

Antibacterial

FOOT CARE PRODUCTS



*Kills up to 99.9% of bacteria and fungus**



*ASTM E-2180-07 Test Method- MRSA, E.coli, Trichophyton rubrum, Trichophyton mentagrophytes

Corn Pads



- ♦ Antibacterial formula
- ♦ Helps to relieve painful corns and calluses
- ♦ Medical grade mineral oil gel pad cushions, protects, and moisturizes
- ♦ Relieves pressure and friction
- ♦ Latex free and hypoallergenic

ITEM #	SIZE	QTY
6513	Medium	6/polybag
6514	Wide	6/polybag

Digital Pads



- ♦ Antibacterial formula
- ♦ Protects toes and fingers from friction and abrasion
- ♦ Helps to alleviate footwear pressure on toes
- ♦ Moisturizes, soothes, comforts, and protects
- ♦ Use where corns, calluses, or hammertoes cause discomfort
- ♦ Latex free and hypoallergenic

ITEM #	SIZE	QTY
6503	Mini	6/polybag
6504	S/M	6/polybag
6505	L/XL	6/polybag


Digital Caps



- ◆ Antibacterial formula
- ◆ Relieves pressure and friction
- ◆ Relieves painful ingrown toenails, corns, hammertoes, and toe tips
- ◆ Fully lined with medical grade mineral oil gel to surround, soothe, moisturize, and protect
- ◆ Fits toes or fingers and can be cut to size
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
6506	S/M	6/polybag
6507	L/XL	6/polybag

Digital Pads on a Strip



- ◆ Antibacterial formula
- ◆ Medical grade mineral oil gel pads on a continuous strip
- ◆ Fabric cuts easily to preferred size
- ◆ After sizing, product should be turned inside out so gel is applied directly to the skin
- ◆ Provides a more cost effective alternative to Digital Pads
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
6500	Mini	17 pads per strip, 2 strips/polybag
6501	S/M	16 pads per strip, 2 strips/polybag
6502	L/XL	14 pads per strip, 2 strips/polybag

Gel Toe Spreaders



- ◆ Antibacterial formula
- ◆ Aligns big toe to help reduce pain
- ◆ Effective on bunions and overlapping toes
- ◆ Separates toes to relieve pressure and friction
- ◆ Used as an inter-digital wedge for bunion deformity, overlapping toes, or toe drift
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
6508	S	6/polybag
6509	M	6/polybag
6510	L	6/polybag

Gel Toe Separators



- ◆ Antibacterial formula
- ◆ Relieves pressure and friction between toes
- ◆ Anatomically designed to fit comfortably between toes
- ◆ Gently divides, cushions, and aligns toes
- ◆ Helps correct overlapping toes
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
6511	S	6/polybag
6512	L	6/polybag

All Gel Bunion Cushions




- ◆ Antibacterial formula
- ◆ Provides relief from bunions
- ◆ Soft, flexible all gel shield cushions, protects, and reduces pressure
- ◆ Thin design fits most footwear
- ◆ Helps to reduce scarring after surgery
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
6515	Hallux Bunion Cushion (Big toe)	1/pkg
6516	Tailor's Bunion Cushion (Little toe)	1/pkg



Corn Pads



- ◆ Helps to relieve painful corns and calluses
- ◆ Medical grade mineral oil gel pad cushions, protects, and moisturizes
- ◆ Relieves pressure and friction
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
1116	Narrow	6/polybag
1117	Medium	6/polybag
1118	Wide	6/polybag

Digital Caps



- ◆ Relieves pressure and friction
- ◆ Relieves painful ingrown toenails, corns, hammertoes, and toe tips
- ◆ Fully lined with medical grade mineral oil gel to surround, soothe, moisturize, and protect
- ◆ Fits toes or fingers and can be cut to size
- ◆ Latex free and hypoallergenic

 **RIBBED
KNIT FABRIC**

 **MESH FABRIC**
(for greater compression)

ITEM #	SIZE	QTY
10195	Mini	6/polybag
10205	S/M	6/polybag
10206	L/XL	6/polybag
10207	XXL	6/polybag

ITEM #	SIZE	QTY
10235	S/M	8/polybag

Digital Pads



- ◆ Protects toes and fingers from friction and abrasion
- ◆ Helps to alleviate footwear pressure on toes
- ◆ Moisturizes, soothes, comforts, and protects
- ◆ Use where corns, calluses, or hammertoes cause discomfort
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
10055	Mini	12/polybag
10065	S/M	12/polybag
10075	L/XL	12/polybag

Digital Pads on a Strip



- ◆ Medical grade mineral oil gel pads on a continuous strip
- ◆ Fabric cuts easily to preferred size
- ◆ After sizing, product should be turned inside out so gel is applied directly to the skin
- ◆ Provides a more cost effective alternative to our Digital Pads
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
1105	Mini	17 pads per strip, 2 strips/polybag
1106	S/M	16 pads per strip, 2 strips/polybag
1107	L/XL	14 pads per strip, 2 strips/polybag

Digistrip



- ◆ Coated on one side with medical grade mineral oil gel strip
- ◆ Fabric cuts easily to preferred size
- ◆ After sizing, product should be turned inside out so gel is applied directly to the skin
- ◆ Provides a more cost effective alternative to our Corn Pads
- ◆ Approximately 30 pads per strip (when cut to 3cm lengths)
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
1100	Narrow, 24" (60cm) long	2 strips/polybag
1101	Medium, 24" (60cm) long	2 strips/polybag
1102	Wide, 24" (60cm) long	2 strips/polybag

Digital Cap with Toe Spreader



- ◆ Traditional Digital Cap with added benefit of built-in Gel Toe Spreader
- ◆ Fully lined with medical grade mineral oil gel for all-around digital protection
- ◆ Relieves pressure and friction
- ◆ Comfortable and simple solution for patients suffering from multiple ailments
- ◆ Aids in the relief of Functional Hallux Limitus when turned with the spreader on the bottom of the toe
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
1596	L/XL	2/polybag
1597	XXL	2/polybag



Gel Tubing



- ◆ Fully lined with medical grade mineral oil gel to provide all around digital protection
- ◆ Less bulky than foam tubing and will not bottom out
- ◆ Can aid in reducing scar tissue
- ◆ Cut to desired size
- ◆ Latex free and hypoallergenic



TRIBBED KNIT FABRIC

ITEM #	SIZE	QTY
10865	Narrow, 3/4" (2cm) diam., 6" (15cm) length	2/polybag
10875	Wide, 1" (2.5cm) diam., 6" (15cm) length	2/polybag



MESH FABRIC (for greater compression)

ITEM #	SIZE	QTY
10855	Narrow, 3/4" (2cm) diam., 6" (15cm) length	3/polybag
10845	Wide, 3" (7.6cm) diam., 10" (25.4cm) length	1/polybag

All Gel Digital Caps



- ◆ Relieves pressure and friction
- ◆ Helps reduce discomfort caused by corns, blisters, and ingrown nails
- ◆ Cushions, soothes, and moisturizes
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
10495	S/M	2/pkg
10505	L/XL	2/pkg

Gel Toe Spreaders



- ◆ Aligns big toe to help reduce pain
- ◆ Effective on bunions and overlapping toes
- ◆ Separates toes to relieve pressure and friction
- ◆ Used as an inter-digital wedge for bunion deformity, overlapping toes, or toe drift
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
11505	Small	4/polybag
11515	Medium	4/polybag
11525	Large	4/polybag

Gel Toe Separators



- ◆ Relieves pressure and friction between toes
- ◆ Anatomically designed to fit comfortably between toes
- ◆ Gently divides, cushions, and aligns toes
- ◆ Helps correct overlapping toes
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
11605	Small	15/polybag
11615	Large	15/polybag

Gel Toe Crest



- ◆ Helps to relieve pressure and pain caused by hammer, claw, and mallet toes
- ◆ Absorbs shock and provides even pressure distribution
- ◆ Reinforced toe loop to ensure proper fit
- ◆ Will not bottom out like foam products
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
10415	Small, Left	3/polybag
10425	Small, Right	3/polybag
10435	Large, Left	3/polybag
10445	Large, Right	3/polybag
10515	Small, Left	1/box
10525	Small, Right	1/box
10535	Large, Left	1/box
10545	Large, Right	1/box

Digital Care Kit



- ◆ A convenient and economical way to store, display, and dispense our most popular digital care products
- ◆ Special compartments for each product allows you to easily select the best solution for the condition you are treating

ITEM #	SIZE
30365	1 kit

KIT INCLUDES

ITEM	SIZE	QTY
Digital Caps	S/M	6
Digital Pads	Mini	4
	S/M	4
	L/XL	4
Mesh Tubing	3/4" x 6"	2
Gel Toe Separators	Small	6
	Large	6
Gel Toe Spreaders	Medium	4
	Large	4

All Gel Care Kit



- ◆ Offers an assortment of our most popular all gel foot care products in a compact kit for easy in-office dispensing
- ◆ An excellent compliment to our Digital Care Kit

ITEM #	SIZE
31744	1 kit

KIT INCLUDES

ITEM	SIZE	QTY
All Gel Digital Caps	S/M	4
	L/XL	4
Gel Metatarsal Pad w/ Toe Spreader	S/M	Pair
	L/XL	Pair
Gel Toe Separators	Small	6
	Large	6
Gel Toe Spreaders	Medium	4
	Large	4



All Gel Hallux/Tailor's Bunion Shield



- ◆ Provides relief from bunions
- ◆ Soft, flexible all gel shield cushions, protects, and reduces pressure
- ◆ Thin design fits most footwear
- ◆ Helps to reduce scarring after surgery
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
1119	Hallux Bunion Shield (Big toe)	1/polybag
1120	Tailor's Bunion Shield (Little toe)	1/polybag

Gel Bunion



- ◆ Anatomically designed to comfortably conform to the shape of the bunion
- ◆ Medical grade mineral oil gel pad cushions, protects, and reduces pressure on affected area
- ◆ Helps to reduce scarring after surgery
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
10085	One Size	5/pkg

Bunion Care Gel Sleeve



- ◆ Anatomically designed to comfortably conform to the shape of the bunion
- ◆ Cushions, protects, and reduces pressure on affected area with medical grade mineral oil gel pad
- ◆ Helps to reduce scarring after surgery
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
10325	S/M	1/polybag
10335	L/XL	1/polybag

Slim Gel-Fit™ Bunion Sleeve



- ◆ Elastic bunion sleeve with the gel pad incorporated inside two layers of a comfortable and conforming fabric
- ◆ The covered gel pad shields the painful and sensitive bony protrusion of the forefoot
- ◆ Helps to relieve shoe pressure, friction, and shear forces on the bunion
- ◆ May also be worn on the hand to protect the thumb from discomfort during gripping activities
- ◆ Thin and compact design makes it perfect for everyday use and fits with most dress, casual, or athletic shoes and underneath hosiery
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
1110	Small	1/polybag
1111	Medium	1/polybag

Universal Gel Strap



- ◆ Comfortable elastic sleeve with anatomically designed gel pad to cover all five metatarsal heads
- ◆ Reduces pain and discomfort in the forefoot
- ◆ Serves as a soft tissue supplement for fat pad atrophy
- ◆ Multi-functional use; cut stitching to fit to instep, bunions, elbow, hand, or wrist
- ◆ Available with uncovered (exposed) gel pad or covered gel pad
- ◆ Latex free and hypoallergenic

GEL UNCOVERED

Mineral oil gel pad helps soften and reduce callused skin

GEL COVERED

Extended product life when callusing is not an issue

ITEM #	SIZE	QTY	ITEM #	SIZE	QTY
10155	S/M	1 pair/polybag	10175	S/M	1 pair/polybag
10165	L/XL	1 pair/polybag	10185	L/XL	1 pair/polybag

Slim Gel-Fit™ Metatarsal Pad



- ◆ Elastic gel metatarsal sleeve with the gel pad incorporated inside two layers of a comfortable and conforming fabric
- ◆ Reduces pain under the ball of the foot caused by metatarsalgia, calluses, rheumatism, neuromas, or from long hours on your feet
- ◆ Thin and compact design makes it perfect for everyday use and fits with most dress, casual, or athletic shoes and underneath hosiery
- ◆ Great for use with high heels that put extra pressure on the metatarsal heads
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
1113	Small	1 pair/polybag
1114	Medium	1 pair/polybag

Gel Foot Cover



- ◆ Ideal for activities that place stress on toes and metatarsals
- ◆ Medical grade mineral oil gel covered with fabric offers complete toe protection
- ◆ Extends to provide protection and cushioning to the metatarsal heads
- ◆ Designed to be worn comfortably with everyday footwear and underneath hosiery
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
10275	Small, 3.75" x 4.9" (9.5cm x 12.5cm)	1 pair/polybag
10285	Medium, 4" x 5.3" (10cm x 13.5cm)	1 pair/polybag
10295	Large, 4.3" x 5.9" (11cm x 15cm)	1 pair/polybag

Ball of Foot Gel Cushion



- ◆ Fabric covered gel pad protects against abrasion, friction, and pressure
- ◆ Ideal for calluses and small neuromas
- ◆ Excellent for relieving the burning sensation under the metatarsal heads
- ◆ May be used as a soft tissue supplement for metatarsalgia
- ◆ Thin design fits most footwear
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
10025	One Size	2/polybag

Gel Metatarsal Pad



- ◆ Absorbs shock, impact, and provides even distribution of pressure in the forefoot area
- ◆ Ideal for conditions such as fat pad atrophy, metatarsalgia, and Morton's neuroma
- ◆ Anatomically designed for comfort and support
- ◆ Helps to reduce the burning sensation in the metatarsal heads
- ◆ Silicone construction with reinforced toe loop to ensure proper fit
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
10465	One Size	1 pair/pkg

Gel Metatarsal Pad with Toe Spreader



- ◆ Helps reduce pain and burning sensations under the metatarsal heads
- ◆ Ideal for fat pad atrophy, metatarsalgia, and Morton's neuroma
- ◆ Anatomically designed for comfort and protection of the metatarsal heads
- ◆ Silicone construction with built-in toe spreader helps reduce discomfort between toes while securing product in place
- ◆ Ideal for activities that place stress on toes and metatarsals
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
4020	S/M	1 pair/box
4021	L/XL	1 pair/box

Sandal Gel Toe Protector



- ◆ Small and sleek gel cushion fits easily around toe stand of sandals and flip flops
- ◆ Cushions between toes and prevents friction for comfortable all day use
- ◆ Peel and stick bottom helps keep gel cushion in place
- ◆ Clear and slim design is barely noticeable
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
1121	One Size	1 pair/polybag

Gel Ball of Foot Pad



- ◆ Absorbs shock and pressure on the ball of foot
- ◆ Gently massages the bottom of the foot and provides extra cushioning
- ◆ Helps prevent calluses by reducing friction and abrasion
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
1582	One Size	1 pair/polybag

Achilles Heel Pad



- ◆ Soft conforming gel pad protects the posterior of the heel from pressure, impact, and friction
- ◆ Provides gentle compression and dissipation of pressure around inflamed area
- ◆ Effective in protecting against discomfort caused by ski boots, skates, cleats, other athletic footwear, shoes, and boots
- ◆ Reduces pain caused by footwear pressure on irregular heel bumps
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
10385	S/M	1/polybag
10395	L/XL	1/polybag

Malleolar Sleeve



- ◆ Cushions and protects the ankle against impact, pressure, and abrasion
- ◆ Designed with two discs of gel molded on both the medial and lateral malleolae
- ◆ Provides gentle compression to the ankle
- ◆ Effective in protecting against discomfort caused by ski boots, skates, cleats, and other athletic footwear
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
15005	One Size	1/polybag

Boot Bumper



- ◆ Helps protect the foot and ankle from abrasions without the added bulk
- ◆ Effective in protecting against discomfort caused by ski boots, skates, cleats, and other athletic footwear
- ◆ Provides light compression to the foot and ankle
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
1770	S/M	1/pkg
1772	L/XL	1/pkg

Skate Bite Protector



- ◆ Anatomically designed gel pad protects the dorsal aspect of the foot from pressure, friction, and pinching from shoelaces
- ◆ Effective in protecting against discomfort caused by ski boots, skates, cleats, and other athletic footwear
- ◆ Provides light compression to the foot and ankle
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
92573	One Size	1/polybag



Soft Skin Heel Sleeve



- ◆ Designed to moisturize and soften dry, cracked skin and helps to reduce heel fissures
- ◆ Thin gel pad anatomically molded to fit the heel
- ◆ Beneficial for dermatitis and eczema
- ◆ Can be used in conjunction with skin moisturizing creams and lotions
- ◆ Recommended as a non-ambulatory therapy product, not to be worn with shoes
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
15265	One Size	1/polybag

Heel/Elbow Slip Over



- ◆ Aids in the healing and skin conditioning of dry, cracked heels and elbows
- ◆ Versatile medical grade mineral oil gel pad can be positioned almost anywhere around the foot and heel
- ◆ Effective in helping to prevent heel and elbow ulcerations
- ◆ Heel sleeve can be used in conjunction with pressure relief products to protect the skin
- ◆ Elbow sleeve provides Ulnar Nerve protection
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
15225	S/M	1/polybag
15235	L/XL	1/polybag
15245	XXL	1/polybag

Half Coated Body Sleeving



- ◆ Moisturization and mild compression attributes make this product excellent for skin grafts, burns, scar treatment, and Lymphedema applications
- ◆ Knit sleeving with moisturizing medical grade mineral oil gel coating on one side of the sleeve
- ◆ 10 foot rolls can be cut to desired length
- ◆ Excellent for use on upper and lower extremities
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
10925	3" x 10' (8.3cm x 3.1m)	1 roll/polybag
10935	4" x 10' (10cm x 3.1m)	1 roll/polybag
10945	5" x 10' (13cm x 3.1m)	1 roll/polybag

Ambidextrous Carpal Gel Sleeve



- ◆ Two thumb holes allow the product to be universally worn on either the left or right hand
- ◆ Medical grade mineral oil gel pad moisturizes and conditions the skin increasing mobility and elasticity of the skin
- ◆ Soft gel pad protects the scarred tissue from friction and pressure
- ◆ Helps to improve and reduce the appearance of scar tissue
- ◆ Reduces discomfort caused by hypersensitivity
- ◆ Product may also be chilled and used for Cold Therapy
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
14165	One Size, 9" to 10" (22.8cm to 25.4cm)	1/polybag

Carpal Gel Sleeve



- ◆ Anatomically designed gel pad molded to elastic sleeve contours to the shape of the palm
- ◆ Helps prevent and reduce shock, impact, and vibration
- ◆ Reduces discomfort caused by hypersensitivity
- ◆ For placement directly over carpal tunnel post-operative scar area
- ◆ Medical grade mineral oil gel pad helps prevent and reduce scar tissue
- ◆ Can be worn under gloves and/or splints
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
14115	Small, Left	1/polybag
14125	Small, Right	1/polybag
14135	Large, Left	1/polybag
14145	Large, Right	1/polybag

Gel-E Roll



- ◆ Elastic bandage with one side coated in a thin layer of medical grade mineral oil gel
- ◆ Conforms to difficult or irregular surfaces of the body and will not slide
- ◆ Gel moves with the skin and does not promote skin irritation
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
10355	4" x 48" (10cm x 120cm)	1/polybag

- May be used:**
- to reduce scar tissue
 - for cold therapy
 - to provide compression
 - as a residual limb shrinker
 - to aid in muscle strains or arthritic pain
 - to moisturize areas of dry skin
 - for burn care
 - to help control edema

Post Wound/Burn Care Kit



- ◆ An assortment of our most popular medical grade mineral oil gel products in a portable kit
- ◆ Indications: Skin breakdown, ulcers, turf toe, blisters, and pressure sores
- ◆ At the first site of pressure sore, place gel over skin to prevent further breakdown
- ◆ Once wound is closed, place gels on skin to help prevent re-opening
- ◆ Gel products diffuse medical grade mineral oil that provides moisturization and elasticity to the skin to promote healing
- ◆ Gel moves independently so it will not cause friction like silicone

ITEM #	QTY
30325	1 kit

KIT INCLUDES

ITEM	SIZE	QTY
10235 - Mesh Digital Cap	S/M	8
10195 - Digital Caps	Mini	6
10205 - Digital Caps	S/M	6
10206 - Digital Caps	L/XL	6
10207 - Digital Caps	XXL	6
10855 - Gel Tubing- Mesh	3/4" x 6"	3
10845 - Gel Tubing- Mesh	3" x 10"	1
10865 - Gel Tubing- Ribbed	3/4" x 6"	2
10875 - Gel Tubing- Ribbed	1" x 6"	2
10305 - Body Discs	2.5"	2
10315 - Body Discs	4"	2
10355 - Gel E Roll	4" x 48"	1
10935 - Half Coated Body Sleeving	4" x 10' roll	1
15235 - Heel/Elbow Slip Over	L/XL	2
Re-Order Card		1

WonderSpur



- ◆ Dual density silicone provides maximum support and shock absorption
- ◆ Aids in calcaneal spurs, plantar fasciitis, heel spurs, and knee & back pain
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	MENS	WOMENS	QTY
4500	Small	2-4	4-5	1 pair/box
4501	Medium	5-7	6-9	1 pair/box
4502	Large	8-10	10-12	1 pair/box
4503	X-Large	11+	12+	1 pair/box

Wondercup



- ◆ Dual density silicone provides relief to entire heel
- ◆ Aids in bursitis, achilles tendonitis, pressure ulcers, and plantar fasciitis
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	MENS	WOMENS	QTY
4504	Small	2-4	4-5	1 pair/box
4505	Medium	5-7	6-9	1 pair/box
4506	Large	8-10	10-12	1 pair/box
4507	X-Large	11+	12+	1 pair/box

WonderSport - Firm



- ◆ Multi-density silicone design simultaneously supports heel and arch
- ◆ Aids in relief for diabetic fat pad atrophy
- ◆ Latex free and hypoallergenic

ITEM #	MENS	QTY	ITEM #	WOMENS	QTY
4508	5-7	1 pair/box	4513	4-5	1 pair/box
4509	8-9	1 pair/box	4514	6-7	1 pair/box
4510	10-11	1 pair/box	4515	8-9	1 pair/box
4511	12-13	1 pair/box	4516	10+	1 pair/box
4512	14+	1 pair/box			

WonderSole Full Length Insole - Firm



- ◆ Full length silicone insole designed to provide support to metatarsals, arch, and heel
- ◆ Helps relieve symptoms related to metatarsalgia, plantar fasciitis, fat pad atrophy, and arthritis
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	MENS	WOMENS	QTY
4530	Small	7	5-7	1 pair/box
4531	Medium	8-9	8-9	1 pair/box
4532	Large	10	10-11	1 pair/box
4533	X-Large	11-13	--	1 pair/box

Heel Spur with Removable Plug - Firm



- ◆ Shock absorbing silicone heel cup with removable plug
- ◆ Wear without plug to allow for pressure disbursement around heel spur
- ◆ Wear with plug for general heel shock absorption
- ◆ Fits all shoe styles
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	MENS	WOMENS	QTY
4001	Small	--	5-6.5	1 pair/box
4002	Medium	5-7.5	7-9.5	1 pair/box
4003	Large	8-11.5	10-12	1 pair/box

Low Profile Heel Cup



- ◆ Soft silicone heel cup provides relief to the entire heel and evenly distributes weight
- ◆ Low profile design fits comfortably in all shoes and provides maximum shock absorption
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	MENS	WOMENS	QTY
4004	Small	--	5-6.5	1 pair/box
4005	Medium	5-7.5	7-9.5	1 pair/box
4006	Large	8-11.5	10-12	1 pair/box

Softzone 3/4 Length Insole - Soft



- ◆ Multi-density silicone anatomically designed to support the longitudinal and transverse arches
- ◆ Deep heel cup ensures superior comfort and control
- ◆ Dual density provides shock absorption and pressure relief to sensitive areas
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	MENS	WOMENS	QTY
4010	Small	--	4.5-6	1 pair/box
4011	Medium	5-7.5	6.5-9.5	1 pair/box
4012	Large	8-10	10-11.5	1 pair/box

Softzone Full Length Insole - Soft



- ◆ Softer multi-density silicone absorbs shock under high pressure areas while supporting the rest of the foot
- ◆ Ideal as a replacement for shoes with removable insoles
- ◆ Relief for heel, arch, and forefoot problems
- ◆ Multiple vents enhance air circulation
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	MENS	WOMENS	QTY
4013	Small	5-6.5	6.5-8.5	1 pair/box
4014	Medium	7-9.5	9-11	1 pair/box
4015	Large	10-12	11.5-12	1 pair/box
4016	X-Large	12.5	--	1 pair/box

Low Profile Fabric-Covered Insole



- ◆ General shock absorption, pressure distribution, and added comfort
- ◆ Soft, durable, and comfortable thin fabric-covered insole allows foot to slide easily into shoe
- ◆ May be used as a lining under other insoles for additional shock absorption and protection
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	MENS	WOMENS	QTY
4030	Small	--	4.5-7	1 pair/box
4031	Medium	5.5-9.5	7.5-11	1 pair/box
4032	Large	10-12.5	11.5-12	1 pair/box

Arthritic/Diabetic Gel Sock



- ◆ Provides optimum protection for the neuropathic foot
- ◆ Non-restrictive cuff
- ◆ Can be used with or without diabetic shoes
- ◆ Helps reduce friction, abrasion, and shear forces while preventing callus formation
- ◆ Single layer sock with 1/16" (2mm) of gel throughout the sole
- ◆ Helps prevent foot ulcerations
- ◆ Serves as a soft tissue supplement for fat pad atrophy
- ◆ Latex free and hypoallergenic

☐ WHITE

ITEM #	SIZE	QTY
1701	US Size 7-9 (UK Size 37-40) Shoe: Womens 6 or smaller	1 pair/box
1702	US Size 9-11 (UK Size 41-43) Shoe: Womens 6 and up Mens 9 or smaller	1 pair/box
1703	US Size 11-13 (UK Size 44-46) Shoe: Mens 9 and up	1 pair/box

☐ BLACK

ITEM #	SIZE	QTY
1706	US Size 7-9 (UK Size 37-40) Shoe: Womens 6 or smaller	1 pair/box
1707	US Size 9-11 (UK Size 41-43) Shoe: Womens 6 and up Mens 9 or smaller	1 pair/box
1708	US Size 11-13 (UK Size 44-46) Shoe: Mens 9 and up	1 pair/box

SoftSock



- ◆ Indicated for neuropathic and rheumatoid arthritis sufferers
- ◆ Excellent for any condition that requires a soft, friction absorbing foot bed
- ◆ Fluid like properties of the gel help massage and comfort the foot while walking
- ◆ Effective in reducing shock and impact on the heel
- ◆ Double layer sock with 1/4" (6mm) of gel throughout the sole
- ◆ Serves as a soft tissue supplement for fat pad atrophy
- ◆ Requires extra depth footwear
- ◆ Latex free and hypoallergenic

☐ WHITE

ITEM #	SIZE	QTY
10605	Small, US Size 7-9 (UK Size 37-40)	1 pair/polybag
10615	Medium, US Size 9-11 (UK Size 41-43)	1 pair/polybag
10625	Large, US Size 10-13 (UK Size 43-46)	1 pair/polybag
10635	X-Large, US Size 13-16 (UK Size 46+)	1 pair/polybag

Gel Dots



- ◆ Self-adhesive backing allows for use with most materials (not skin)
- ◆ Gel dot is 1/8" (3mm) thick and 1" (3cm) in diameter
- ◆ Latex free and hypoallergenic

ITEM #	QTY
10255	15 dots/sheet/polybag

- May be used:**
- to reduce abrasion, friction, and pressure points in splints, casts, braces, or footwear
 - can be fabricated with insoles and placed under individual metatarsal heads
 - to protect the ulnar styloid
 - as individual cushions for parts of the hand such as the knuckles in casts or splints

Gel Squares



- ◆ Versatile solution for discomfort related to pressure and friction
- ◆ Gel square is 1/8" (3mm) thick and 4" x 4" (10cm x 10cm) in size
- ◆ Can be trimmed to customized size
- ◆ Medical grade mineral oil gel moisturizes and conditions the skin
- ◆ Can be used as an interface between the skin and a brace
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
15505	With Self-Adhesive Backing	2 pieces/polybag
15515	Without Backing	2 pieces/polybag

Body Discs



- ◆ Multi-purpose discs of medical grade mineral oil gel
- ◆ Apply to parts of the body suffering from discomfort caused by pressure and friction
- ◆ Can be trimmed to customized size
- ◆ Medical grade mineral oil gel moisturizes and conditions the skin
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
10305	2.5" (6.3cm) Diameter	2 pieces/polybag
10315	4" (10cm) Diameter	2 pieces/polybag

Carpal Gel Pads



- ◆ Anatomically designed gel pad (same as in the Carpal Gel Sleeve) with self-adhesive backing
- ◆ Can be incorporated into custom made or pre-fabricated hand and wrist braces
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
14205	2" x 6" (5cm x 15.24cm)	3 pieces/polybag

Pressure Relief Padding



- ◆ Used to reduce shock, impact, and vibration
- ◆ Can be used as a replacement for moleskin or felt products
- ◆ Excellent for use with custom ankle foot orthosis devices
- ◆ Perfect for cushion protection in hand splints
- ◆ 1/16" (2mm) thick gel pad with self-adhesive backing
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
11205	4" x 36" (10cm x 91cm)	1/pkg
11215	8" x 36" (20cm x 91cm)	1/pkg

Gel and Hook Brace Pads



- ◆ 1/4" (6mm) gel thickness to protect skin from friction and shearing
- ◆ Attaches to foam on inside of brace
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
1785	Circle, 2.5" (6.3cm) Diameter	2/polybag
31740	Rectangle, 1.75" x 2.5" (4.5cm x 6.3cm)	1/polybag

Soft Shear GEL Sheeting



- ◆ 1/8" (3mm) thick gel sheet available with medical grade mineral oil gel exposed on one side, or lined with fabric on both sides
- ◆ Fabric top cover provides excellent stretch properties
- ◆ Relieves pressure and shear forces
- ◆ Fabric side can be glued to braces or other materials
- ◆ Latex free and hypoallergenic


FABRIC ONE SIDE

ITEM #	SIZE	QTY
15605	8" x 8" (20cm x 20cm)	1/pkg
15705	12" x 16" (30.5cm x 40.5cm)	1/pkg
15805	16" x 20" (40.5cm x 51cm)	1/pkg

FABRIC TWO SIDES

ITEM #	SIZE	QTY
15615	8" x 8" (20cm x 20cm)	1/pkg
15715	12" x 16" (30.5cm x 40.5cm)	1/pkg
15815	16" x 20" (40.5cm x 51cm)	1/pkg

WonderFlex SILICONE Sheets



- ◆ Silicone sheeting available in multiple thicknesses
- ◆ Available lined or unlined
- ◆ Wonderflex sheeting is 12" x 16" (30.5cm x 40.5cm) in size
- ◆ Latex free and hypoallergenic

LINED

ITEM #	SIZE	QTY
4733	3mm	1/pkg
4735	5mm	1/pkg

- Soft durometer cushioning for top cover of orthotics
- Ideal for diabetic, arthritic, and sensitive feet
- Resists compression

UNLINED

ITEM #	SIZE	QTY
4722	2mm	1/pkg
4723	3mm	1/pkg

- Recommended for use with orthotic devices



Cushion, protect, and relieve foot pain and pressure with every step. TheraStep® foot care products turn uncomfortable walks into soothing therapy sessions for troubled feet.

Each product relieves the pain associated with bunions, corns, blisters, and calluses comfortably and effectively, as medical grade mineral oil moisturizes and softens the skin.

The TheraStep® consumer packaged line offers our best selling products for professional in-office dispensing. Increase revenue for your practice while providing your patients with a proven solution for their ailments.

TheraStep® products are latex free and hypoallergenic.

Gel Corn Pads



ITEM #	SIZE	QTY
6205	One Size	2/pkg

Gel Toe Cap



ITEM #	SIZE	QTY
6203	Small	1/pkg
6204	Large	1/pkg

Gel Toe Tubes



ITEM #	SIZE	QTY
6208	One Size	2/pkg

All-Gel Toe Cap



ITEM #	SIZE	QTY
6201	Small	1/pkg
6202	Large	1/pkg

Gel Toe Spreaders



ITEM #	SIZE	QTY
6207	Medium	2/pkg

Gel Toe Separators



ITEM #	SIZE	QTY
6206	2 Small 2 Large	4/pkg

Cushion Toe Separators



ITEM #	SIZE	QTY
6216	2 Small 2 Medium 2 Large	6/pkg

Cushion Toe Comb



ITEM #	SIZE	QTY
6220	One Size	2/pkg

Gel Toe Support



ITEM #	SIZE	QTY
6212	Left, One Size	1/pkg
6213	Right, One Size	1/pkg

Gel Bunion Shield



ITEM #	SIZE	QTY
6209	One Size	1/pkg

Cushion Ball of Foot Pad



ITEM #	SIZE	QTY
6215	One Size	2/pkg

Gel Ball of Foot Pad



ITEM #	SIZE	QTY
6214	One Size	2/pkg

Heel Cups



ITEM #	SIZE	QTY
6210	One Size	2/pkg

Gel Forefoot Pad



ITEM #	SIZE	QTY
6211	One Size	2/pkg

Arch Support



ITEM #	SIZE	QTY
6224	One Size	2/pkg





Prosthetics





Designed with the patient in mind, Silipos' prosthetic liners and suspension sleeves provide users with an unmatched level of comfort and confidence. Unlike silicone, our medical grade mineral oil gels move with the skin and will not pull or irritate the residual limb. Moisturizing medical grade mineral oils are released from our gels, aiding in the reduction of scar tissue and helping with overall skin conditioning. We offer a wide array of products that are sure to meet the needs of your patients.

Our prosthetic products are field tested with industry experts in order to assure patient comfort and satisfaction.

Silosheath™

Silosheath™ provides optimal range of motion and comfort, and can be used with other interfaces as necessary. Made with Silipos medical grade mineral oil gel and nylon outer material, Silosheath™ is available in 3 different styles:

SILOSHEATH™ SIZING GUIDELINES

SIZE	DESCRIPTION	GEL LENGTH	LENGTH	DISTAL CIRCUMFERENCE	PROXIMAL CIRCUMFERENCE
					
1	Child	8" (20cm)	10"-13" (25cm-33cm)	6.5"-9" (17cm-23cm)	6.5"-10" (17cm-25cm)
2	Short/Narrow	10" (25cm)	15"-17" (38cm-43cm)	7"-10.5" (18cm-27cm)	7"-12" (18cm-30.5cm)
3	Short/Regular	10" (25cm)	15"-17" (38cm-43cm)	8"-12.5" (20cm-32cm)	8"-14" (20cm-35.5cm)
4	Medium/Narrow	13" (33cm)	17"-19" (43cm-48cm)	8"-12.5" (20cm-32cm)	8"-14" (20cm-35.5cm)
5	Medium/Regular	10" (25cm)	17"-19" (43cm-48cm)	9"-14" (23cm-35.5cm)	9"-16" (23cm-40.5cm)
6	Large/Narrow	13" (33cm)	18"-21" (46cm-53cm)	9"-14" (23cm-35.5cm)	9"-16" (23cm-40.5cm)
7	Large/Regular	13" (33cm)	18"-21" (46cm-53cm)	11"-16" (28cm-41cm)	11"-18" (28cm-46cm)
8	Large/Wide	13" (33cm)	18"-21" (46cm-53cm)	12"-19" (30.5cm-48cm)	12"-21" (30.5cm-53cm)

All sizes are approximate.

All Silosheath™ products offer the following benefits:

Durability

Conformability

Maximum Protection

Unmatched patient comfort and satisfaction

Can be used as a night time skin conditioner and post op scarring aid

Latex free and hypoallergenic

Silosheath™ Original



◆ 2 ply, 2mm gel

ITEM #	SIZE	QTY
13315	1 (Child)	1/box
13325	2 (Short/Narrow)	1/box
13335	3 (Short/Regular)	1/box
13345	4 (Medium/Narrow)	1/box
13355	5 (Medium/Regular)	1/box
13365	6 (Large/Narrow)	1/box
13375	7 (Large/Regular)	1/box
13385	8 (Large/Wide)	1/box

Silosheath™ Double Cushion



◆ 2 ply, 3mm gel, Extra gel at distal end

ITEM #	SIZE	QTY
13415	1 (Child)	1/box
13425	2 (Short/Narrow)	1/box
13435	3 (Short/Regular)	1/box
13445	4 (Medium/Narrow)	1/box
13455	5 (Medium/Regular)	1/box
13465	6 (Large/Narrow)	1/box
13475	7 (Large/Regular)	1/box
13485	8 (Large/Wide)	1/box

Silosheath™ Extra Life







◆ 3 ply, 2mm gel, Increased durability

ITEM #	SIZE	QTY
13815	1 (Child)	1/box
13825	2 (Short/Narrow)	1/box
13835	3 (Short/Regular)	1/box
13845	4 (Medium/Narrow)	1/box
13855	5 (Medium/Regular)	1/box
13865	6 (Large/Narrow)	1/box
13875	7 (Large/Regular)	1/box
13885	8 (Large/Wide)	1/box

ComfortZone™ Gel Liners

Silipos Gel Liners offer the highest level in quality, protection, and comfort. ComfortZone™ liners offer maximum protection with 3mm thick medical grade mineral oil gel. Available in 3 different styles:

COMFORTZONE™ SINGLE SOCKET GEL LINER SIZING GUIDELINES

ITEM#	SIZE	PLY	GEL LENGTH	LENGTH OF LINER	DISTAL CIRCUMFERENCE	PROXIMAL CIRCUMFERENCE
						
12725	Short/Narrow	3	10" (25cm)	12" (30.5cm)	9"-12" (23cm-30.5cm)	15"-18" (38cm-45cm)
12755	Short/Regular	3	10" (25cm)	12" (30.5cm)	12"-14" (30.5cm-35.5cm)	15"-18" (38cm-45cm)
12765	Medium/Regular	3	12" (30.5cm)	16" (40.5cm)	12"-14" (30.5cm-35.5cm)	15"-18" (38cm-45cm)

All sizes are approximate.

COMFORTZONE™ DOUBLE SOCKET GEL LINER SIZING GUIDELINES

ITEM#	SIZE	PLY	GEL LENGTH	LENGTH OF LINER	DISTAL CIRCUMFERENCE	PROXIMAL CIRCUMFERENCE
12625	Short/Narrow	5	10" (25cm)	12" (30.5cm)	9"-12" (23cm-30.5cm)	15"-18" (38cm-45cm)
12655	Short/Regular	5	10" (25cm)	12" (30.5cm)	12"-14" (30.5cm-35.5cm)	15"-18" (38cm-45cm)
12665	Medium/Regular	5	12" (30.5cm)	16" (40.5cm)	12"-14" (30.5cm-35.5cm)	15"-18" (38cm-45cm)

All sizes are approximate.

COMFORTZONE™ ULTRA CUSHION GEL LINER SIZING GUIDELINES

ITEM#	SIZE	PLY	GEL LENGTH	LENGTH OF LINER	DISTAL CIRCUMFERENCE	PROXIMAL CIRCUMFERENCE
12555	I-Short/Regular	3	10" (25cm)	12.5" (32cm)	7"-11" (18cm-28cm)	8"-13" (20cm-33cm)
12515	I-Short/Regular	5	10" (25cm)	12.5" (32cm)	7"-11" (18cm-28cm)	8"-13" (20cm-33cm)
12565	II-Medium/Regular	3	12" (30.5cm)	15"-17" (38cm-43cm)	8"-12.5" (20cm-32cm)	8"-14" (20cm-35.5cm)
12525	II-Medium/Regular	5	12" (30.5cm)	17"-19" (43cm-48cm)	8"-12.5" (20cm-32cm)	8"-14" (20cm-35.5cm)
12575	III-Medium/Wide	3	12" (30.5cm)	17"-19" (43cm-48cm)	9"-14" (23cm-35.5cm)	9"-16" (23cm-40.5cm)
12535	III-Medium/Wide	5	12" (30.5cm)	18"-21" (46cm-53cm)	9"-14" (23cm-35.5cm)	9"-16" (23cm-40.5cm)

All sizes are approximate.

All ComfortZone™ Gel Liner products offer the following benefits:

- Prevents abrasion and skin breakdown
- Optimum cushion and protection
- Ease of mobility
- Unmatched patient comfort and satisfaction
- Latex free and hypoallergenic

Single Socket Gel Liner



- ◆ 3 ply, 3mm of exposed gel, CoolMax® exterior fabric, tapered shape

ITEM #	SIZE	PLY	QTY
12725	Short / Narrow	3	1/box
12755	Short / Regular	3	1/box
12765	Medium / Regular	3	1/box

Double Socket Gel Liner



- ◆ 5 ply, 3mm gel, two layers of CoolMax® exterior fabric, tapered shape

ITEM #	SIZE	PLY	QTY
12625	Short / Narrow	5	1/box
12655	Short / Regular	5	1/box
12665	Medium / Regular	5	1/box

Ultra Cushion Gel Liner



- ◆ 3 ply or 5 ply, 3mm gel, two layers of cotton material, nylon-reinforced distal end

ITEM #	SIZE	PLY	QTY
12555	I - Short / Regular	3	1/box
12565	II - Medium / Regular	3	1/box
12575	III - Medium / Wide	3	1/box
12515	I - Short / Regular	5	1/box
12525	II - Medium / Regular	5	1/box
12535	III - Medium / Wide	5	1/box



SiloLiner™ is the most skin-friendly and easy to use solution for geriatric and post-op patients.

- ◆ Easy to Donn and Doff: Reinforced cuff prevents liner from rolling down and irritating the skin
- ◆ Medical grade mineral oil gel helps moisturize and condition skin grafts and scars
- ◆ Outer material provides wide range of motion
- ◆ 6mm gel thickness
- ◆ Latex free and hypoallergenic

SILOLINER™ PIN LINER SIZING GUIDELINES

ITEM#		FLAT WIDTH	DISTAL CIRCUMFERENCE at urethane attachment	PROXIMAL CIRCUMFERENCE MPT
1460	1464	3" (7.6cm)	8"-12" (20cm-30.5cm)	8"-14" (20cm-35.5cm)
1461	1465	4" (10.2cm)	9.5"-13.5" (24cm-34cm)	9.5"-16" (24cm-40.5cm)
1462	1466	5" (12.7cm)	12"-15" (30.5cm-38cm)	12"-20" (30.5cm-51cm)
1463	1467	6" (15.2cm)	14"-16" (35.5cm-40.5cm)	14"-26" (35.5cm-66cm)

All sizes are approximate.

SILOLINER™ CUSHION LINER SIZING GUIDELINES

ITEM#		FLAT WIDTH	DISTAL CIRCUMFERENCE at urethane attachment	PROXIMAL CIRCUMFERENCE MPT
1612	1622	3" (7.6cm)	8"-14" (20cm-35.5cm)	8"-14" (20cm-35.5cm)
1613	1623	4" (10.2cm)	9.5"-16" (24cm-40.5cm)	9.5"-16" (24cm-40.5cm)
1614	1624	5" (12.7cm)	12"-17" (30.5cm-43cm)	12"-20" (30.5cm-51cm)
1615	1625	6" (15.2cm)	14"-19" (35.5cm-48cm)	14"-26" (35.5cm-66cm)

All sizes are approximate.

SiloLiner™ Pin Liner



PIN LINER

ITEM #	LENGTH	FLAT WIDTH	QTY
1460	12" (30.5cm)	3" (7.6cm)	1/box
1461	12" (30.5cm)	4" (10.2cm)	1/box
1462	12" (30.5cm)	5" (12.7cm)	1/box
1463	12" (30.5cm)	6" (15.2cm)	1/box
1464	14" (35.5cm)	3" (7.6cm)	1/box
1465	14" (35.5cm)	4" (10.2cm)	1/box
1466	14" (35.5cm)	5" (12.7cm)	1/box
1467	14" (35.5cm)	6" (15.2cm)	1/box

SiloLiner™ Cushion



CUSHION LINER

ITEM #	LENGTH	FLAT WIDTH	QTY
1612	12" (30.5cm)	3" (7.6cm)	1/box
1613	12" (30.5cm)	4" (10.2cm)	1/box
1614	12" (30.5cm)	5" (12.7cm)	1/box
1615	12" (30.5cm)	6" (15.2cm)	1/box
1622	14" (35.5cm)	3" (7.6cm)	1/box
1623	14" (35.5cm)	4" (10.2cm)	1/box
1624	14" (35.5cm)	5" (12.7cm)	1/box
1625	14" (35.5cm)	6" (15.2cm)	1/box

Low Activity Liner

The Low Activity Liner provides superior skin conditioning and comfort for low to moderate activity level patients.

- ♦ Stabilizing Matrix reduces distal elongation and minimizes pistoning
- ♦ Liner easily conforms to limbs
- ♦ Retrofittable to Silipos Explorer™ sockets
- ♦ Available in 3mm, 6mm, and 9mm gel thickness
- ♦ Latex free and hypoallergenic

Low Activity Pin Liner



PIN LINER, 3MM GEL THICKNESS

ITEM #	SIZE	STYLE	QTY
18101	Medium	Uniform	1/box
18107	Medium+	Uniform	1/box
18113	Large	Uniform	1/box
18116	Large+	Uniform	1/box
18119	XLarge	Uniform	1/box

ITEM #	SIZE	STYLE	QTY
18102	Medium	Tapered	1/box
18108	Medium+	Tapered	1/box

PIN LINER, 6MM GEL THICKNESS

ITEM #	SIZE	STYLE	QTY
18103	Medium	Uniform	1/box
18109	Medium+	Uniform	1/box
18114	Large	Uniform	1/box
18117	Large+	Uniform	1/box
18120	XLarge	Uniform	1/box



ITEM #	SIZE	STYLE	QTY
18104	Medium	Tapered	1/box
18110	Medium+	Tapered	1/box

PIN LINER, 9MM GEL THICKNESS

ITEM #	SIZE	STYLE	QTY
18105	Medium	Uniform	1/box
18111	Medium+	Uniform	1/box
18115	Large	Uniform	1/box
18118	Large+	Uniform	1/box

ITEM #	SIZE	STYLE	QTY
18106	Medium	Tapered	1/box
18112	Medium+	Tapered	1/box

LOW ACTIVITY LINER SIZING GUIDELINES

PROXIMAL CIRCUMFERENCE			DISTAL CIRCUMFERENCE	
 Measure 15" from Distal End			 Measure 1.5" from Distal End	
SIZE	RELAXED	FULLY STRETCHED	RELAXED	FULLY STRETCHED
Medium	9" (23cm)	16" (41cm)	8" (20cm)	12" (30cm)
Medium+	13" (33cm)	21" (53cm)	9" (23cm)	12" (30cm)
Large	15" (38cm)	22" (56cm)	9" (23cm)	15" (38cm)
Large+	16" (41cm)	27" (68cm)	12" (30cm)	16" (41cm)
X-Large	16" (41cm)	30" (76cm)	12" (30cm)	20" (51cm)

All sizes are approximate.

Low Activity Cushion Liner



CUSHION LINER, 3MM GEL THICKNESS

ITEM #	SIZE	STYLE	QTY
18141	Medium	Uniform	1/box
18147	Medium+	Uniform	1/box
18153	Large	Uniform	1/box
18156	Large+	Uniform	1/box
18159	XLarge	Uniform	1/box

ITEM #	SIZE	STYLE	QTY
18142	Medium	Tapered	1/box
18148	Medium+	Tapered	1/box

CUSHION LINER, 6MM GEL THICKNESS

ITEM #	SIZE	STYLE	QTY
18143	Medium	Uniform	1/box
18149	Medium+	Uniform	1/box
18154	Large	Uniform	1/box
18157	Large+	Uniform	1/box
18160	XLarge	Uniform	1/box

ITEM #	SIZE	STYLE	QTY
18144	Medium	Tapered	1/box
18150	Medium+	Tapered	1/box

CUSHION LINER, 9MM GEL THICKNESS

ITEM #	SIZE	STYLE	QTY
18145	Medium	Uniform	1/box
18151	Medium+	Uniform	1/box
18155	Large	Uniform	1/box
18158	Large+	Uniform	1/box

ITEM #	SIZE	STYLE	QTY
18146	Medium	Tapered	1/box
18152	Medium+	Tapered	1/box



DuraGel™ offers an unparalleled combination of durability, comfort, and skin protection for moderate to above moderate activity level patients.

- ◆ Nylon composite fabric provides extreme abrasion resistance and superior gel bond
- ◆ Gel formulation will not tear, delaminate, or thin out with stress, maintaining skin care and overall patient comfort
- ◆ Retrofittable to Silipos Explorer™ sockets
- ◆ Available in 3mm, 6mm, and 9mm gel thickness
- ◆ Latex free and hypoallergenic

DuraGel™ Pin Liner



PIN LINER, 3MM GEL THICKNESS

ITEM #	SIZE	STYLE	QTY
18211	Medium	Uniform	1/box
18212	Medium+	Uniform	1/box
18213	Large	Uniform	1/box
18214	Large+	Uniform	1/box
18215	XLarge	Uniform	1/box

ITEM #	SIZE	STYLE	QTY
18216	Medium	Tapered	1/box
18217	Medium+	Tapered	1/box

PIN LINER, 6MM GEL THICKNESS



ITEM #	SIZE	STYLE	QTY
18201	Medium	Uniform	1/box
18202	Medium+	Uniform	1/box
18203	Large	Uniform	1/box
18204	Large+	Uniform	1/box
18205	XLarge	Uniform	1/box

ITEM #	SIZE	STYLE	QTY
18206	Medium	Tapered	1/box
18207	Medium+	Tapered	1/box
18208	Large	Tapered	1/box

PIN LINER, 9MM GEL THICKNESS

ITEM #	SIZE	STYLE	QTY
18221	Medium	Uniform	1/box
18222	Medium+	Uniform	1/box
18223	Large	Uniform	1/box
18224	Large+	Uniform	1/box

DURAGEL™ LINER SIZING GUIDELINES

PROXIMAL CIRCUMFERENCE			DISTAL CIRCUMFERENCE	
 Measure 15" from Distal End			 Measure 1.5" from Distal End	
SIZE	RELAXED	FULLY STRETCHED	RELAXED	FULLY STRETCHED
Medium	9" (23cm)	15" (38cm)	8" (20cm)	13" (33cm)
Medium+	13" (33cm)	20" (51cm)	9" (23cm)	13" (33cm)
Large	15" (38cm)	21" (53cm)	9" (23cm)	16" (41cm)
Large+	16" (41cm)	26" (66cm)	12" (30cm)	17" (43cm)
X-Large	16" (41cm)	29" (74cm)	12" (30cm)	21" (53cm)

All sizes are approximate.

DuraGel™ Cushion Liner



CUSHION LINER, 3MM GEL THICKNESS

ITEM #	SIZE	STYLE	QTY
18251	Medium	Uniform	1/box
18252	Medium+	Uniform	1/box
18253	Large	Uniform	1/box
18254	Large+	Uniform	1/box
18255	XLarge	Uniform	1/box

ITEM #	SIZE	STYLE	QTY
18256	Medium	Tapered	1/box
18257	Medium+	Tapered	1/box

CUSHION LINER, 6MM GEL THICKNESS


ITEM #	SIZE	STYLE	QTY
18241	Medium	Uniform	1/box
18242	Medium+	Uniform	1/box
18243	Large	Uniform	1/box
18244	Large+	Uniform	1/box
18245	XLarge	Uniform	1/box

ITEM #	SIZE	STYLE	QTY
18246	Medium	Tapered	1/box
18247	Medium+	Tapered	1/box
18248	Large	Tapered	1/box

CUSHION LINER, 9MM GEL THICKNESS

ITEM #	SIZE	STYLE	QTY
18261	Medium	Uniform	1/box
18262	Medium+	Uniform	1/box
18263	Large	Uniform	1/box
18264	Large+	Uniform	1/box

DuraGel™ Suspension Sleeve




The DuraGel™ Suspension Sleeve combines the same gel used in our DuraGel™ liners with a flexible outer fabric to provide improved durability and range of motion.

- ◆ Skin friendly gel works well with suction application
- ◆ Easy to Donn and Doff
- ◆ DuraGel™ mineral oil gel moisturizes and conditions the skin while providing increased durability
- ◆ 3mm gel thickness
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
18181	Medium	1/polybag
18182	Medium+	1/polybag
18185	Medium+ Distal+	1/polybag
18183	Large	1/polybag
18184	Large+	1/polybag

DuraGel™ Reinforced Suspension Sleeve



Internal fabric reinforcement provides increased durability without limiting flexibility. This sleeve maintains suspension and provides added security for its wearer.

ITEM #	SIZE	QTY
18191	Medium	1/polybag
18192	Medium+	1/polybag
18195	Medium+ Distal+	1/polybag
18193	Large	1/polybag
18194	Large+	1/polybag

DURAGEL™ SUSPENSION SLEEVE SIZING GUIDELINES

SIZE	PROXIMAL CIRCUMFERENCE		DISTAL CIRCUMFERENCE	
	RELAXED	FULLY STRETCHED	RELAXED	FULLY STRETCHED
Medium	9" (23cm)	15" (38cm)	8" (20cm)	13" (33cm)
Medium+	14.25" (36cm)	26.75" (68cm)	9.25" (23cm)	16" (40cm)
Medium+ Distal+	14.25" (36cm)	26.75" (68cm)	9.75" (24.5cm)	17" (43cm)
Large	13.5" (34cm)	24" (61cm)	12" (31cm)	18" (46cm)
Large+	16.25" (41cm)	28" (71cm)	11.25" (29cm)	20" (51cm)

All sizes are approximate.

BK Suspension Sleeve



- ◆ Provides maximum suspension
- ◆ Lightweight and durable
- ◆ Ease of application
- ◆ Seamless construction
- ◆ Medical grade mineral oil gel moisturizes and conditions the skin
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	APPROX. LENGTH	APPROX. CIRCUMFERENCE	QTY
13055	I	12" (30.5cm)	11"-17" (28cm-43cm)	1/polybag
13065	II	14" (35.5cm)	13"-20" (32.5cm-51cm)	1/polybag
13075	III	16" (40.5cm)	14"-22" (35.5cm-56cm)	1/polybag
13085	IV	16" (40.5cm)	15"-28" (38cm-71cm)	1/polybag

HD-BK Suspension Sleeve



- ◆ Helps to prevent rotation
- ◆ Internal fabric reinforcement increases durability
- ◆ Will not bunch-up behind the popliteal area
- ◆ Ease of application
- ◆ Medical grade mineral oil gel moisturizes and conditions the skin
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	APPROX. LENGTH	APPROX. CIRCUMFERENCE	QTY
1310	I	12" (30.5cm)	11"-17" (28cm-43cm)	1/polybag
1311	II	14" (35.5cm)	13"-20" (32.5cm-51cm)	1/polybag
1312	III	16" (40.5cm)	14"-22" (35.5cm-56cm)	1/polybag

Arm Gel Suspension Sleeve



- ◆ Compatible with myo-electric and cosmetic prosthesis
- ◆ May also be used as BK sleeve for children and small adults
- ◆ Medical grade mineral oil gel moisturizes and conditions the skin
- ◆ Latex free and hypoallergenic

ITEM #	APPROX. LENGTH	APPROX. CIRCUMFERENCE	QTY
13005	8"-13" (20cm-30cm)	8.75" (22cm)	1/polybag

Symes Gel Liner



- ◆ 5 ply, 3mm gel
- ◆ Retrofittable to most Symes' socket designs
- ◆ Prevents abrasion and skin breakdown
- ◆ Optimum conformability
- ◆ Made with medical grade mineral oil gel
- ◆ CoolMax® Material provides elasticity, comfort, and breathability
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
1302	One Size Fits Most	1/polybag

LENGTH	GEL LENGTH	DISTAL CIRCUMFERENCE	PROXIMAL CIRCUMFERENCE
16.5" (42cm)	15" (38cm)	10"-14.5" (25cm-37cm)	15"-20" (38cm-51cm)

All sizes are approximate.

Partial Foot Sock



- ◆ For use with transmetatarsal amputees
- ◆ CoolMax® Material provides elasticity, comfort, and breathability
- ◆ May be used with foot orthoses or shoes with fillers
- ◆ Medical grade mineral oil gel covers site of amputation
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
10805	Regular	1/polybag
10815	Extra Long	1/polybag

Mineral Oil Gel Roll-On



- ◆ Apply to distal end of limb for added comfort and protection
- ◆ Helps to enhance socket liner fit
- ◆ Adheres to Silipos gel liners to ease application
- ◆ Can be cut to size
- ◆ Medical grade mineral oil gel moisturizes and conditions the skin
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	GEL THICKNESS	APPROX. OUTER CIRC.	APPROX. GEL LENGTH	QTY
1495	Narrow	2mm	7.25" (18.5cm)	6" (15cm)	2/polybag
1496	Medium	2mm	8.5" (21.5cm)	6" (15cm)	2/polybag
1497	Wide	2mm	9.25" (23.5cm)	6" (15cm)	2/polybag

Distal End Pad



- ◆ Provides extra protection for shock and abrasion at the distal end
- ◆ Aids in enhancing socket fit for limbs undergoing changes in volume or mass
- ◆ Can be cut to size
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
12005	One Size, 4" (10cm) Diameter	1/polybag





Skin Care

Silipos gels release moisture-rich medical grade mineral oil that gently hydrates, softens, and conditions the skin. Our gloves and socks are beneficial for treating skin conditions such as eczema, psoriasis, and dermatitis while giving the patient a pampered feeling. Our gel forms protect the post-surgical area while releasing medical grade mineral oils that help to reduce scarring while moisturizing the skin.

Our Gel-Care® Advanced product line offers a revolutionary new scar management system using the clinically proven SymRepair® Technology. The gel sheeting can be cut to desired size and the body wraps offer flexibility for hard to reach, bendable areas of the body. Gel-Care® Advanced is an excellent alternative to silicone for the prevention and reduction of scars.

Silipos gels are washable and reusable, and with proper care will continue to release skin-beneficial mineral oils throughout the life of the product.




Advanced Scar Treatment Products

Gel-Care® Advanced offers a revolutionary new scar management system that delivers advanced ingredients contained in an enriched sheet. The formula continues to treat the skin even after the sheet is removed.

- ◆ Prevents and diminishes the appearance of scars from injury, surgery, burns, acne, and stretch marks
- ◆ Excellent for use in the prevention and treatment of hypertrophic and keloid scars
- ◆ Can be trimmed to customized size
- ◆ Washable and reusable
- ◆ Silicone free, latex free, and hypoallergenic

SymRepair®

Synergistic blend of three key components of the barrier lipid combined with Bisabolol. SymRepair® quickly reconstitutes lamellar bilayers resembling those of the skin.



Coolact 38 D (p-methane 3, 8-diol)

Patented Dermal Promotion Absorption Agent.

High Linoleic Acid (Safflower Oil)

Emollient lipid soothes the skin and helps restore the skin's barrier function.

Retinyl Palmitate (Vitamin A)

Regulates epidermal development and modulates collagen synthesis.

Tocopherol (Vitamin E) Acetate

Antioxidant. Long established as beneficial to proper skin care.

Gel-Care® Advanced Self-Adhesive Sheet




ITEM #	SIZE	QTY
610	5" x 6" (12cm x 15cm)	1/pkg
611	5" x 6" (12cm x 15cm)	10/pkg

Gel-Care® Advanced Self-Adhesive Strips



ITEM #	SIZE	QTY
613	1.5" x 10" (4cm x 25cm)	2/pkg

Gel-Care® Advanced Self-Adhesive Discs



ITEM #	SIZE	QTY
615	1.5" (4cm)	6/pkg

Gel-Care® Advanced Body Wraps



- ◆ Flexibility of application for hard to reach and bendable areas of the body

ITEM #	SIZE	QTY
620	3" x 36" (7.5cm x 91cm)	1/pkg
621	6" x 36" (15cm x 91cm)	1/pkg

Gel-Care®



- ◆ 2mm thick layer of self-adhesive medical grade mineral oil gel with conforming material backing
- ◆ Moisturizes the skin and helps to reduce scarring
- ◆ Washable and reusable
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
760	4" x 4" (10cm x 10cm)	1/box
761	5.5" x 8" (14cm x 20cm)	1/box
762	1.5" x 10" (4cm x 25cm)	2/box

Post Surgical Breast Forms



- ◆ 0.5mm thick layer of self-adhesive medical grade mineral oil gel with conforming nylon backing
- ◆ For use after breast augmentation procedures
- ◆ Allows patients to wear clothing next to hypersensitive skin areas
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
740	Small, 32-36 (81.2cm x 91.4cm)	1 pair/polybag
741	Large, 36-42 (91.4cm x 106.7cm)	1 pair/polybag

Post Surgical Abdominal Form



- ◆ 0.5mm thick layer of self-adhesive medical grade mineral oil gel with conforming nylon backing
- ◆ For use after abdominoplasty procedures
- ◆ Allows patients to wear clothing next to hypersensitive skin areas
- ◆ Latex free and hypoallergenic

ITEM #	SIZE	QTY
744	Small, 17.3" x 2.4" (44cm x 6cm)	1/polybag
745	Large, 17.7" x 3.5" (45cm x 9cm)	1/polybag

Gel Therapy Gloves & Socks



- ◆ Medical grade mineral oil gel comfortably conforms to the hands or feet for continuous moisturization and conditioning of the skin
- ◆ Beneficial for skin conditions such as eczema, psoriasis, and dermatitis
- ◆ Designed with a non-restrictive cuff
- ◆ Can be used as warming gloves/socks for arthritic conditions
- ◆ Gloves and socks can be worn at night for intensive treatment (products are non-ambulatory)
- ◆ Latex free and hypoallergenic and does not contain any additives or fragrances

GEL THERAPY GLOVES

ITEM #	SIZE	COLOR	QTY
786	One Size	White	1 pair/pkg

GEL THERAPY SOCKS

ITEM #	SIZE	COLOR	QTY
788	One Size	Oatmeal	1 pair/box





GeLuscious® features a fun, funky, functional line of well-being and lifestyle products for busy everyday living. Indulge in a sumptuous range of patented products for a luxurious spa experience at home, in the office, or while traveling.

Infused with exclusive vitamin enriched antioxidants, nourishing vitamins, and moisture-rich oils, GeLuscious® polymer gels have been proven to completely transform the skin, making even the driest, most callused skin smooth, supple, and beautiful.

Our footcare gel cushions are a comfortable solution for fashionable but painful shoes.

GL for Men™ gel socks are designed with many pores in their structure to absorb odor and control temperature, while moisturizing gel gets rid of dry, callused skin.

GeLuscious® provides the ultimate experience in comfort, health, and well-being.



GeLuscious® Moisturizing Gel Gloves



ITEM #	COLOR	QTY
271	Turquoise	1 pair/pkg
272	Pink	1 pair/pkg

GeLuscious® Moisturizing Gel Booties



ITEM #	COLOR	QTY
251	Turquoise	1 pair/pkg
253	Pink	1 pair/pkg

GeLuscious® GL for Men™



ITEM #	SIZE	QTY
261	Charcoal	1 pair/pkg

GeLuscious® Smooth Heelin™



ITEM #	QTY
254	1 pair/pkg

GeLuscious® Sassy Strips™



ITEM #	QTY
293	2 pairs/pkg

GeLuscious® Divine Drops™



ITEM #	QTY
294	6/pkg

GeLuscious® Pumpfort Pad™



ITEM #	QTY
291	1 pair/pkg

GeLuscious® Heelmate™



ITEM #	QTY
292	1 pair/pkg

GeLuscious® In-2-Sole™



ITEM #	QTY
295	1 pair/pkg

GeLuscious® Cuti-Fix™



ITEM #	QTY
265	10/pkg and 1 cuticle pen





Market Leading Expertise in Gel Manufacturing, Processing, and Conversion

FDA Registered Medical Device Manufacturing Facility

ISO 9001:2008 and ISO 13485:2003 Certified Company

USP Medical Grade Mineral Oil

Superior Product Innovation and Quality

Made in the USA

Silipos, Inc.
7049 Williams Road
Niagara Falls, NY 14304
United States
Toll Free: 800-229-4404
Phone: +1 (716) 283-0700
Fax: +1 (716) 283-0600

**European Authorized
Representative**
WMDE
Bergerweg 18
6085 AT Horn
Holland

Silipos Asia Pacific
P.O. Box 5122
Mordialloc
Victoria 3195
Australia
Phone: +61 (3) 9590-9818
Fax: +61 (3) 9590-9003

www.silipos.com